



# Sample Preparation Considerations for Laser MicroDissection



# What are the Considerations for Sample Preparation for LCM / LMD?



# Leica LMD Protocol Guide



Living up to Life



## Sample Preparation For Laser Microdissection

Protocol Guide for Leica Microsystems Laser Microdissection Systems

# Staining Protocols



- Less is More
- Adjust protocols for minimal stain to allow identification of target cells.
- Less stain can improve contrast in the image.
- Use no stain.
- Use stained sections on one slide as a navigational tool to identify areas of interest on unstained sections on another slide
- Use software module to transfer the regions of interest from a properly stained section on one slide to unstained serial sections on other slides.
- Fluorescent tags require special consideration.

# Typical Sample Preparation for Laser Microdissection for Frozen Sections for Animal Cells



- Cut sections 4-20 um
- Mount on chosen slides
- Fixation with 70% EtOH for 30 sec. minimum
- Wash in DEPC water for 30 sec.
- Stain with Mayer's hematoxlin for 30 sec.
- Wash in DEPC water for 30 sec.
- Counterstain with eosin for 30 sec.
- Rinse with (95%) and (100%) EtOH for 30 sec.
- Air dry approximately 5 minutes

# Fixation

Preserve the tissues; preserve the macromolecules; stop metabolism; allow visualization identification; better LMD.

## Recommended

- Cryofixation, a good choice
- Fixation by ice-cold acetone, 30 – 45 seconds
- Fixation by graded alcohols up to 70%


## Other

- Xylene fixation for short times. (Xylene substitute not recommended.)
- Formalin fixation @ 4% will cross link RNA,
- Good fixation will allow for better LCM/LMD

# Fluorescence Considerations

- Over-dehydration by Alcohols, or Xylene may cause both immunofluorescence and fluorescent proteins (e.g. GFP) to deactivate.
- Recommend series of alcohol fixation/dehydration up to 70% Ethanol.
- Recommend short time fixation/dehydration in Xylene.
- GFP is a common fluorescent protein tag; it varies widely by manufacturer. Contact with polymers may cause GFP to deactivate. Use a coating (e.g. poly-L-lysine) as a protectant.
- Choose a low auto-fluorescent slide (e.g. PPS, FLC, EP).

# Choose the Right Slide



- Systems using foil slides:
- PET (polyethylene terephthalate) recommended for MALDI, SELDI ionization analysis.
- PEN (polyethylene naphthalate) not recommended for any ionization analysis.
- Expression Pathology EP Director slides recommended for proteomic analysis, for large areas; low auto fluorescence background.
- Low fluorescence background slides with PPS, FLC, EP.



# Pre-treating Slides

- UV irradiation prior to mounting specimen has beneficial effects:
  - Makes all slides RNase, DNase, Protease free and sterile.
  - Makes foil slides more adherent for tissue.
  - Makes laser cutting of foil slides more efficient.

# Know your LCM / LMD System

- What are the capabilities and limitations of the laser for cutting ROI's?
- What is the maximum sample thickness your system is capable of cutting/capturing?
- How can you optimize the system parameters for cutting/capturing your particular specimen?
- How can you save the laser parameters?

# Some References

## Citations:

**Proteomic Analysis of Cells Isolated by Laser Microdissection**, V. Fiorilli, V. Klink, et. al., Chapter 3, “**Integrative Proteomics**”, Intech, 2012, pp 47 – 69.

**Desmin expression in colorectal cancer stroma correlates with advanced stage disease and marks angiogenic microvessels**, G. Arentz, et. al., *Clinical Proteomics* 2011, 8:16

**Increased Phospho-Keratin 8 Isoforms in Colorectal Tumors Associated with EGFR Pathway Activation and Reduced Apoptosis**, G. Arentz, et. al., **Molecular Biology** Volume 2012, Article ID 706545

**MALDI-TOF and MALDI-FTICR imaging mass spectrometry of methamphetamine incorporated into hair**, AkihiroMiki, et. al., *J. Mass. Spectrom.* 2011, 46, 411–416.

**Medullary amyloidosis associated with apolipoprotein A-IV deposition**, S. Sethi, et. al., *Kidney International* 2011 International Society of Nephrology,

# References:

## Websites:

<http://www.leica-microsystems.com/products/light-microscopes/life-science-research/laser-microdissection/details/product/leica-lmd7000/>

[http://www.expressionpathology.com/pdf/iTRAQ\\_labeling\\_of\\_Liquid\\_Tissue\\_Preps\\_Rev1.pdf](http://www.expressionpathology.com/pdf/iTRAQ_labeling_of_Liquid_Tissue_Preps_Rev1.pdf)

<http://www.expressionpathology.com/pdf/LiquidTissueManualV146.pdf>

<http://www.qiagen.com/literature/handbooks/literature.aspx?id=1000939>

<http://www.invitrogen.com/site/us/en/home.html>