Scalable and Highly Sensitive Single Cell RNA-seq through Combinatorial Barcoding

Genomics

Charlie Roco (charlie@parsebiosciences.com), Parse Biosciences, Vuong Tran, Parse Biosciences, Emily OConnor, Parse Biosciences, Peter Matulich, Parse Biosciences, Ryan Koehler, Parse Biosciences, Daniel Diaz, Parse Biosciences, Joey Pangallo, Parse Biosciences, Sarah Schroeder, Parse Biosciences, Grace Wang, Parse Biosciences, Archita Gadkari, Parse Biosciences, Brittany Takushi, Parse Biosciences, Lili Bahrami, Parse Biosciences, Nam Tran, Parse Biosciences, Anna Malinkevich, Parse Biosciences, Moe Blackman, Parse Biosciences, Alex Rosenberg, Parse Biosciences

Combinatorial barcoding has been demonstrated as a promising approach to scaling up single cell and single nuclei transcriptomics. Here, we present a commercial kit that makes it possible to profile up to 48 different biological samples across 100,000 cells in a single experiment. The kit is based on the Split Pool Ligation-based Transcriptome sequencing (SPLiT-seq) protocol, but includes a streamlined workflow and dramatic improvements in performance. Transcript and gene detection have improved over 5-fold compared to the published method. The rates of doublets in resulting data are low even when profiling 100,000 cells (~3%). The workflow also enables biological samples that are collected on different days to be rapidly fixed, stored, and then later pooled and barcoded in parallel in a single experiment. We will share example data collected from various sample types including cell lines, brain, and PBMCs samples.