

## Automated extraction of cfDNA

### Genomics

**Janina Gerhards** (janina.gerhards@perkinelmer.com), PerkinElmer chemagen, **Anne-Maria Wallraf**, PerkinElmer chemagen, **Mareike Honickel**, PerkinElmer chemagen, **Tim Weissbach**, PerkinElmer chemagen

The chemagic cfDNA kits simplify the efficient cfDNA extraction from 0.5-5 ml plasma by enabling an automated workflow on the chemagic 360 instrument with minimal hands-on time. Thereby, the kits show equivalent yields and purities compared to a manual gold standard kit.